

## Establishing New Cultures from Frozen Stocks

### Reagents:

*Note: Sterilize all solutions by 0.22 $\mu$ m filtration and store at 4.0°C.*

- Solution A: 12% w/v NaCl in ddH<sub>2</sub>O
- Solution B: 1.6% w/v NaCl in ddH<sub>2</sub>O
- Solution C: RPMI complete medium with 10% human serum

### Procedure:

1. Remove cryovial containing parasitized red blood cells from cryotank and place at room temperature to defrost. (Cryovial should contain ~500-1000 $\mu$ L)
2. For each 500 $\mu$ L thawed blood, add 100 $\mu$ L of Sol A drop by drop and mix slowly.
3. Transfer the culture to a sterile 15mL tube and let it stand for 3 min. (Move quickly to Step 4)
4. For each original 500 $\mu$ L of thawed blood, add 5mL Sol B drop by drop and mix slowly by pipetting up and down.
5. Centrifuge at 2000 rpm (800xg), 5min at RT. (Don't forget to set the brake to low)
6. Remove and discard supernatant.
7. For each original 500 $\mu$ L of thawed blood, add 5mL Sol C.
8. Centrifuge at 2000rpm (800xg), 5min at RT.
9. Remove and discard the supernatant.
10. Resuspend the culture in 10 mL of Sol C in a small culture flask, adding 500 $\mu$ L of fresh washed 50% hematocrit blood to obtain a final hematocrit of ~5%.
11. Gas the culture (1% O<sub>2</sub>, 3%CO<sub>2</sub>, 96% N<sub>2</sub>) and place in 37°C incubator.
12. Change the medium and smear daily to monitor the parasitemia. When the parasitemia reaches 8-10%, the stock culture should be diluted down and large cultures may then be set up using complete medium with 5% human serum and Albumax.